


DENTAL LINES

DENTAL UNITS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Luxury Dental Chair

Model Series YR06012



Product Description

Kalstein YR06012 dental comprehensive treatment, which is upgraded comprehensively, is easy to operate, has excellent quality, complete functions and exquisite appearance, which gives you double enjoyment of vision and touch. Perfect patient chair and doctor chair structure design, bring patients and doctors enjoy the experience. The new reconstruction of the arm space structure brings luxury diagnosis and treatment space. Free all the troubles for you.

- Enjoy unprecedented comfortable experience. The structural design and manufacturing process are in line with the most comfortable ergonomics and the strictest safety standards. The independent patient seat is designed with the best proportion of linkage to ensure that patients enjoy the most comfortable experience.
- Luxuriant diagnosis space for treat. It is the best choice for high-class clinic and VIP Clinic because of its wide range mobile arm system, easy to achieve the ergonomic positioning of dentists, fashionable appearance and luxurious treatment space design.
- Safety and health. The back of chair and ceramic spittoon stop when they are blocked, and the handpiece and chair interlock function can ensure the safety of the whole treatment process.
- Dental led light with 8 lights, there are cold light and warm light, touch screen, can also be manually controlled color temperature between 4000-5500K, illumination up to 60000 Lux.
- Reserved loading position of microscope.
- Hygienic ceramic spittoon with stainless steel handle, flushing and water supply pipe are convenient for cleaning and disinfection.
- Flexible option. There are three models for you to choose. Hanging , top mounted and cart model. For cart model, The height adjustment range is 650-850mm, with brake silent wheel. The top mounted type provides a wide treatment range for dentist, and the maximum extension distance of the instrument can reach 88mm.
- Exquisite structure and fashionable appearance meet the aesthetic requirements of dental clinic.
- Expand the operation field of vision and control unit, improve work efficiency.
- The multifunctional integrated foot switch controls all chair positions and spittoons ,free dentist hands.
- The super skin doctor's chair, ergonomic chair surface and four-way adjustable back support can relieve your tension and fatigue.
- The chair mounted microscope can be optioned.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Luxury Dental Chair

Model Series YR06012



Standard Configuration:

Model	YR06012
Leather cushion seat	1 set
3-way syringe(cold/hot)	2 sets
Powerful suction device	1 set
saliva ejector	1 set
Dentists Stool (Leather cushion)	1 set
Warm water system	1 set
Pure water supply system	1 set
LED x-ray viewer	1 set
Rotatable and removeable spittoon	1 set
DC 24V motor driven chair	1 set
Multifunctional Assistant control	1 set
Synchronized chair movement	1 set
Dental unit Foot pad	1 set
Multifunctional foot control pedal	1 set
Backrest, spittoon anti-collision protection function	1 set
Luxury 8 LED Sensor light	1 set



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Luxury Dental Chair

Model Series YR06012



Technical Description::

Model	YR06012
Power voltage	220V-50Hz /110V-60Hz
Motor voltage	24V
Air pressure	0.5mpa-0.8mpa
Water pressure	0.2mpa-0.4mpa



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most used equipment in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog


Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

