



LABORATORY SECTOR

MICROSCOPE



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Microscope
Biological Microscope
Model YR0242-M // YR0242-T



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Microscope Biological Microscope Model YR0242-M // YR0242-T



Product Description:

- ✓ Monocular, binocular, trinocular and teaching head available for choose.
- ✓ Quadruple nosepiece with precise click stop.
- ✓ Integrated double layer mechanical stage, size 142x132mm, moving range 75x50mm.
- ✓ Abbe condenser N.A.1.25, high brightness 3w led, brightness adjustable
- ✓ Infinity plan and achromatic optical system available

Seidentopf Head:

Monocular, binocular, trinocular and teaching head available for choose. All heads inclined at 30°, rotatable 360°. For binocular and trinocular, both eyepiece tube diopter adjustable +/-5°.

Quadruple, Backward:

Quadruple nosepiece with precise click stop. Equipped with 35mm achromatic CF objective 4X, 10X, 40Xs, 100Xs (oil).

Double Slide Holder Working Stage:

Integrated double layer mechanical stage, size 142x132 mm, moving range 75x50 mm.

Condenser and Light Source:

Abbe condenser N.A.1.25 high brightness 3W LED, brightness adjustable.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Microscope Biological Microscope Model YR0242-M // YR0242-T



Technical Specifications:

Model No.	YR0242-M	YR0242-B	YR0242-T
Video Adapter	—	C-Mount 0.3x	—
Power Supply Type	AC 110V-240V	—	—
Transmit Light	Halogen 6V20W or LED Optional		
Light Source Type	Halogen		
Condenser N.A.	1.25		
Condenser Type	Abbe		
Condenser Features	Rack & Pinion↓↑, Filter Holder, Center Adjustable		
Focusing Scale	0.002mm		
Focusing Range	36mm		
Focusing Type	Coaxial Coarse & Fine	—	—
Focusing Features	Safety Stop Screw Tension Adjustable		
Working Stage Moving	75*50mm		
Working Stage Size	132*142mm		
Working Stage Type	2 Layers Mechanical		
Objective	Achromatic 4x, 10x, 40x(S), 100x(O,S)		
Objective Type	Achromatic		
Nosepiece Holes	4 Holes		
Nosepiece Features	Inward		
Eyepiece	WF10x/18mm		
Eyepiece Diameter	23.2mm		
Head Interpupillary	54~75mm		
Head Rotation	Head Rotate 360°	—	—
Head Inclination	30°	—	—
Head Type	Monocular	Binocular	Seidentopf Trinocular
Standard Magnification	40x-1000x	—	—
Optical System	Achromatic		
Standard Function	BF	—	—



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

