




LABORATORY LINES

COLOR LIGHT BOX



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



CONTROL THE EVER-CHANGING COLORS

When testing product color or color matching in production process, different light sources, illuminance, color temperature and environmental changes make the naked eye make wrong judgment on color. Color light box can provide accurate and objective standard light source environment for raw material production or product color comparison.



PRODUCT DESCRIPTION



F



D65



TL84



UV



CWF



TL83



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Color Light Box

Model YR05520 // YR05523



Technical specifications

Model	YR05520	YR05521	YR05522	YR05523	
Light source	D65, TL84, UV, F, CWF, TL83	D65, TL84, UV, F, CWF	D65, TL84, UV, F, CWF, TL83	D65, TL84, UV, F, CWF	
Dimensions(W*D*H)	695x552x502 mm	695x552x502 mm	695x552x502 mm	695x552x502 mm	
Inner Size(W*D*H)	665x490x360 mm				
Voltage	220V/50Hz – 110V/60Hz				
Light sources Configuration	D65	2 PCS	2 PCS	2 PCS	2 PCS
	TL84	2 PCS	2 PCS	2 PCS	2 PCS
	F	4 PCS	4 PCS	4 PCS	4 PCS
	UV	1 PCS	1 PCS	1 PCS	1 PCS
	CWF	2 PCS	2 PCS	2 PCS	2 PCS
	TL83	2 PCS	/	2 PCS	/
Weight	19.5 kg	19.5 kg	28.9 kg	28.9 kg	
Light Source Specifications					
Light Source	Color Light Box Application	Color temperature(k)	Power(W)		
D65	International standard Artificial Daylight	6500	18		
TL84	Applied to stores in Europe, Japan and China	4000	18		
CWF	(Cool White Fluorescent)American Standard	4150	18		
F	“Sun-setting Light Yellow” incandescent light source (imitation of sunset)	2700	40		
UV	Viewing under ultraviolet light to detect and evaluate optical brighteners or fluorescent pigments	Wavelength:365nm	20		
TL83	Applied to stores in Europe	3000	18		
Remarks: Optional tubes D50, D75, U30, U35, A					



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



COLOR LIGHT BOX



Model YR05520

1. Composite engineering plastics, die production, no yellowing, deformation, paint removal;
2. Seamless structural design, can add elastic pad, light baffle (optional);
3. Displays the display time, name and switching times of each light source;
4. The light source can be switched freely, with the function of same color different spectrum;
5. No preheating, low energy consumption, no heating, high luminous efficiency;
6. More complete British and American standard common light source, light source name can be changed;
7. Size: 69.5*55.2*50.2cm;
8. Light source: D65 TL84 CWF UV F TL83.



Model YR05526 Color Light Box (Adjustable light lux)
Light source : F, D65, TL84, UV, CWF, TL83



Model YR05529 Color light box
Light Source : D65, TL84, TL83, UV, F, CWF six light sources



Model YR05536 Color Light Table (Color proof table)
Light Source Optional: D50(5000K) or D65(6500K) or TL84 (4000K);

CC120 Series Color Light Table (Color proof table)



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



COLOR LIGHT BOX



Model YR05530D

Material: HDF

Features: Equipped with common standard light source, suitable for multiple light source environment;

Light Source configuration: A,F, D65, TL84, UV ,CWF,TL83/U30 four light sources;

Color temperature type: International brand lamp, standard gray environment;

Size: 710*530*570mm

Application: Color management in textile, toys, printing and dyeing, plastics, paint, ink, printing, pigment, chemical industry, ceramics, shoes, leather, hardware, food, cosmetics and other industries

MORE MODELS AND CUSTOM LIGHT SOURCES

Model	Light source	D65	D50	TL84	TL83/U30	UV	F	CWF	A
YR05520	Six light source	●		●	TL83	●	●	●	
YR05521	Five light source	●		●		●	●	●	
YR05522	Six light source	●		●	TL83	●	●	●	
YR05523	Five light source	●		●		●	●	●	
YR05526	Six light source	●		●	TL83	●	●	●	
YR05524	Four light source	●		●		●	●		
YR05525	Five light source	●		●		●	●	●	
YR05527	Four light source	●		●		●	●		
YR05528	Six light source	●		●	U30	●		●	●
YR05529	Six light source	●		●	TL83	●	●	●	
YR05530	Seven light source	●		●	TL83	●	●	●	●
YR05530	Eight light source	●		●	TL83/U30	●	●	●	●
YR05530A	Four light source	●		●		●	●		
YR05530B	Five light source	●		●		●	●	●	
YR05530C	Six light source	●		●	TL83	●	●	●	
YR05530D	Seven light source	●		●	TL83	●	●	●	●
YR05533		●	●	●					
YR05534		●	●	●					



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



STANDARD LIGHT SOURCE

Provide the following Lamps:

D65, D50, D75, TL84, CWF, UV, U30, TL83, U35, F, A, INCA, HOR



3nhLighting D65 lamp
 Model: F20T12/65 6500K 20W
 Size: 60cm
 Brand: 3nh



PHILIPS D65 lamp
 Model: TLD18W/965
 Size: 60cm
 Brand: PHILIPS



GRETAGMACBETH D65 lamp
 Model: F20T12/65 6500K
 Size: 60cm
 Brand: Macbeth



VeriVide D65 lamp
 Model: F20T12/D65
 Size: 60cm
 Brand: VeriVide



Ecolux U35 lamp
 Model: F17T8 SPX35 ECO
 Size: 60cm
 Brand: GE



SYLVANIA D65 lamp
 Model: F20T12/65 6500K 20W
 Size: 60cm
 Brand: SYLVANIA

Name	Color temperature	Kalstein	PHILIPS	VeriVide	GE	SYLVANIA	GretagMacbeth
D75 North Sky Daylight	7500K	✓					✓
D65 Artificial Daylight	6500K	✓	✓	✓		✓	✓
D50 Professional light source for printing	5000K	✓	✓				✓
CWF Cool White Fluorescent	4150K		✓		✓	✓	
TL84 Commercial lighting applications outside North America	4000K	✓		✓			
HOR Horizontal fluorescent lamp	3200K	✓				✓	
TL83 Warm White Fluorescent	3000K		✓			✓	
A American window spotlights, colorimetric reference lights	2856K	✓			✓	✓	
F Family hotel lights, colorimetric reference light source	2700K	✓					
UV Ultra-Violet	365nm	✓	✓			✓	



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell, • 34000 Montpellier, France.
 Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Réseaux Sociaux et Sites Officiels)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris – France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu
<https://kalstein.net>
<https://kalstein.eu>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most important pieces of equipment for laboratories. It is a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?


Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Vidéos de formation Kalstein
(Chaîne Youtube et sites officiels)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

