




LABORATORY SECTOR

REFRIGERATORS AND FREEZERS











OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.







FEATURERS



Multi-Scene applications





Heated glass design





Coil evaporator

Better performance, easy maintenance



Suction circulation system
Less temp. flutuation & better



Lockable doors
Standard locks & padlocks



Large LED displays

Easy to observe and operate, built-in USB data



Precise temp. control

Temp. control accurate to 0.1°C

Adjustable shelves and drawers



6 types of alarm

Audible and visual alarms





Variable space







CFC-free refrigerant and energy-saving











PARAMETER

Performance		Refrigeration System		Alarms	
Cooling method	Direct cooling	Refrigerant	R:R600 &F:R290	High/low temperature	Υ
Defrost	Manual	Temp. controller	Microprocessor	Power failure	Υ
Refrigerant	R600a &R290	Sensor	NTC	Sensor error	Υ
Temperature range	F:-20°C~-40°C R:2°C~8°C	Display	Digital display	Controller failure	Υ
Ambient temperature	10~32℃	External material	SPCC sprayed steel	Door ajar	Υ
Noise (db)	50	Condenser	Copper tube	Ambient temperature	Υ
Electrical Data		Dimensions		Other	

Electrical Data		Dime	Dimensions		Other	
Voltage frequency(V/Hz)	220V 50Hz/60Hz	Package size(W*D*H)(mm)	775*920*1990	Inside door	3 shelves & 6 drawers	
Power (W)	430W	Capacity (L/cu.ft)	R:183/6.5 F:185/6.5	Remote alarm port	Optional	
Current (A)	2.2A	Exterior size(W*D*H)(mm)	720*830*1870	Access port QTY/Dia.	2/25MM	
Power consumption (kWh/24h)	6.85	Interior size(W*D*H)(mm)	R:580*626*625 F:466*636*625	Internal material	SPCC sprayed steel (SUS 304 for option)	





#LETSGIVEMORE

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



biomédicale

https://www.humatem.org/





ANY DOUBT?Contact us!

Articles Kalstein

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official website





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier

75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel**: +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/



(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

