



LABORATORY LINES

SHAKER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Decoloration Shaker

Series YR05765



DS Decolorization shaker adopts digital display, dc motor, microcomputer control speed, with power recovery and fault self-diagnosis.

DS Decolorization shaker is widely used for experiments where samples need a shaking in aqueous solution, fixed electrophoresis gel separation band, shaking while Coomassie blue staining and decolouring, silver nitrate staining fixed, dyeing and enhancement.

Advantage

- ✓ Using uniaxial drive balancing technology, machine runs smoothly with low noise.
- ✓ Brushless DC motor, long life, no need maintenance.
- ✓ Overcurrent protection function, automatically cut off power without affecting other equipment.
- ✓ Automatic memory function, after starting up, run according to the last setting
- ✓ Slow start function to prevent liquid splashing

Model	YR05765	YR05765-1	YR05765-2
Shake model	convolution/rectiprocation	convolution/rectiprocation	convolution/rectiprocation
Cyclotron frequency range	40-300rpm	40-300rpm	40-300rpm
Frequency sensitivity(rpm)	±1	±1	±1
Orbit Diameter(mm)	Φ26	Φ26	Φ26
Standard capacity	9x250ml	12x250ml	25x250ml
Size of platform(mm)	330x280	350x350	490x490
Timing range (continuous available)	0-990min	0-990min	0-990min
Number of platform	1	1	1
Outside dimension (W*D*H) (mm)	350x330x230	350x400x230	510x540x230
Power (W)	300	400	500
Power supply	AC220V±10% 50-60HZ	AC220V±10% 50-60HZ	AC220V±10% 50-60HZ
Optional	Double layer, max capacity will be double of standard capacity		
	Max capacity for conical/erlenmeyer flask		
50ml	20	36	49
100ml	16	28	39
250ml	9	12	25
300ml	8	12	20
500ml	6	10	16
1L	3	6	9
2L	/	4	6



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

