




LABORATORY SECTOR

INCUBATORS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Incubators
Dry Bath Incubator
Model YR05895 // YR05896



#Letsgivingmore ♥


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Incubators Dry Bath Incubator Model YR05895 // YR05896



Product Description:

YR05894 is a dry block heater controlled by microcomputer. It can be widely used in sample preservation and reaction, DNA amplification and pre-denaturation of electrophoresis, serum coagulation, etc. in fields of medicine, chemical industry, food safety, quality control and environment.

Features:

- ✓ LED display with timer
- ✓ High precision temperature
- ✓ Built-in over-temperature protect
- ✓ Small size with transparent lid
- ✓ Various blocks can protect samples from contamination



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Incubators Dry Bath Incubator Model YR05895 // YR05896



Technical Specifications:

Model	YR05895	YR05896
Temperature control mode	Semi-conductor	Heating membrane
Temperature control range	RT-30°C ~100°C	RT+5°C ~100°C
Temperature accuracy	±0.3°C	
Temperature Stability @40°C	±0.3°C	
Temperature Stability @100°C	±0.3°C	
Timer	0~99h 59m	
Display accuracy	0.1°C	
Heating Up Rate	≤ 15min(20°C to 100°C)	
Max.Temperature	100°C	
Cooling Down Rate	≤ 30min 100-25°C ≤ 30min Cooling down 30degrees from room temperature (test under RT26°C)	—
Power	130W	150W

Optional Accessories:

Type	Capacity	Max. Temperature	Description
YR05896-op1	96×0.2ml PCR plate	100°C	Centrifuge Tube
YR05896-op2	54×0.5ml tubes	100°C	Centrifuge Tube
YR05896-op3	35×1.5ml tubes	100°C	Centrifuge Tube
YR05896-op4	35×2.0ml tubes	100°C	Centrifuge Tube
YR05896-op5	20×0.5ml+15×1.5ml tubes	100°C	Centrifuge Tube
YR05896-op6	24×φ12mm tubes	100°C	Centrifuge Tube
YR05896-op7	32×0.2ml+25×1.5ml tubes	100°C	Centrifuge Tube
YR05896-op8	32×0.2ml+10×0.5ml+15×1.5ml tubes	100°C	Centrifuge Tube
YR05896-op9	103×67×30 (water tank)	100°C	Water Tank
YR05896-op10	96×0.2ml Elisa (Flat bottom)	100°C	Plate (Flat Bottom)
YR05896-op11	24×5ml tubes	100°C	(Speed≤600rpm)
YR05896-op12	12×15ml tubes	100°C	(Speed≤600rpm)
YR05896-op13	6×50ml tubes	100°C	(Speed≤600rpm)
YR05896-op14	12×10ml tubes	100°C	(Speed≤600rpm)
YR05896-op15	Lid	100°C	Used with YR05896-op1–YR05896-op9
YR05896-op16	Lid	100°C	Used with YR05896-op10–YR05896-op11
YR05896-op17	Lid	100°C	Used with YR05896-op12–YR05896-op14



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

