


Kalstein A different accompaniment, at your service

LABORATORY LINES FREEZER DRYER





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Electric Heating Pilot Freezer Dryer

Model Series YR010FA // YR050FD



Electric Heating Pilot Freeze Dryer

Pilot freeze dryer is in situ freeze dryer, it is microprocessor-controlled, pre-cooling and freeze-drying are completed in one process, which is the best choice for scientific research, development and small batch production.

Features

- In Situ Freeze Dryer, Easy to Operate and Good Drying Effect
- Small shelf temperature difference, High temperature control precision and uniform drying effect
- · 7 inch Color Touch Screen, USB interface
- Powerful sensor calibration to ensure accurate measurement of long-term use
- Preserve hundreds of process recipes, each of which contains 36 temperature control sections
- · Flexible manual mode and automatic control mode
- Real-time alarm display and historical alarm guery function
- The user level and password can be set, and decentralized operation management









#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



Electric Heating Pilot Freezer Dryer







Specification

Model	YR010FA	YR010FB	YR030FD	YR050FD
Freeze-dried area	0.2m²	0.1m²	0.4m²	0.6m²
Panel material (material thickness 10mm)	2L	1L	4L	6L
Material tray	395 × 265mm 2 layers	395×265mm 1 layers	295×335mm 4 layers	295×335mm 6 layers
Material plate spacing	70mm			50mm
Shelf temperature range	-50°C~+70°C			
Cold trap temperature	≤ -60°C (no load), optional ≤ -80°C (no load)			
Water catching capacity	≥6KG/24h		≥8KG/24h	≥10KG/24h
Ultimate vacuum	≤ 5pa (no load)			
Power	2000w		2700W	4500w
Host weight	210kg		280kg	450kg
Mainframe dimensions	655×640×1200mm	655×640×1200mm	890×660×1210 mm	890×670×1470mm



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



ANY DOUBT? Contact us!

Articles Kalstein

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

