




MEDICAL SECTOR

WHEELCHAIR



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Wheelchair
Electric Wheelchair
Model YR05445



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Wheelchair Electric Wheelchair Model YR05445



Product Description:

- ✓ Can board
- ✓ Lightweight folding design, easy to carry and travel
- ✓ Intelligent brake
- ✓ Flip-up armrest
- ✓ 360° Adjustable push handle

About US:

Kalstein medical is an integrated industry and trade company focusing on medical device R&D, production, sales and consulting services. At present, the main products that are developed and sold are wheelchairs and medical beds, old-age scooters, walkers, toilet chairs, crutches, etc. We can provide OEM or ODM services according to the different needs of customers. At the same time, we can also recommend products that are popular among local market people based on our global market data analysis and sales. Our factory can provide you with a full range of purely customized services including brand, function, material, parts, etc., to ensure that it can meet the various customized needs of users; Our products are exported to Europe, America, South America, the Middle East and Southeast Asia and other global markets, and we have established in-depth cooperative relations with local distributors, charities, and medical rehabilitation institutions. Our professional technology and industry experience in the medical field are the Professional guarantee for cooperation and service.

Wheelchair:

Kalstein medical wheelchair is our most advantageous product. We can provide wheelchairs with different functions for people with different needs, including electric wheelchairs, manual wheelchairs, wheelchairs for children with cerebral palsy, standing wheelchairs, stair climbing wheelchairs and other wheelchairs with different functions.

Electric Medical Bed:

Kalstein medical electric medical bed is not only able to provide hospital beds, but also home care beds according to personal care needs. In terms of function, there are multiple functions such as electric turning, back lifting, leg lifting, lifting, and defecation to provide options to meet the needs of different scenes and users.



#Letsgivingmore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Wheelchair Electric Wheelchair Model YR05445



Manual Medical Bed

Kalstein medical manual medical bed is a cost-effective product series. According to the needs of different people, there are many functions such as turning over, lifting back, lifting legs, and releasing bowels. At the same time, to meet the needs of hospital beds, we also launched a series of nursing beds for families and individuals, giving customers more choices.

Product Production Process:

- ✓ Purchase and inspection of raw materials and standard parts.
- ✓ Processing and inspection of outsourcing components.
- ✓ Mechanical processing and inspection of parts.
- ✓ Component assembly and inspection.
- ✓ Surface treatment and inspection of parts.
- ✓ Product assembly.
- ✓ Product inspection.
- ✓ Product packaging and product storage.



Note: The surface treatment baking varnish, spraying, and electroplating are all made by outsourcing.

Model	YR05445
Seat width:	36 cm
Seat depth:	36cm
Seat height:	53cm
Backrest height:	43 cm
Overall Length:	68 cm
Overall width:	58cm
Overall height:	35cm
Rear wheel size:	8 inches
Front wheel size:	6 inches
Folding Height:	34 cm
Folding Length:	68 cm
Folding Width:	58 cm
N/G:	18KG/16.8KG
Frame material:	aluminum alloy
Load:	100KG
Driving Range:	12KM
Max Speed:	8KM/h
Motor:	300W*2 Double hub motor
Battery:	24V 6Ah Lithium
Parking Brake:	manual brake and electromagnetic brake
Cushion:	Black



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument used for observing objects too small to be seen by the naked eye. This is achieved through the use of lenses or a mirror and a system of reflecting and refracting the light and other electromagnetic waves from the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device used to store food and other items at a low temperature. It is one of the most common pieces of equipment used in laboratories to maintain a controlled environment (refrigerated space), without fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

