




# LABORATORY SECTOR

## CLIMATIC CHAMBER



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.




### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Climatic Chamber**  
**Environmental Chamber**  
**Model YR05340 // YR05346-1**



**#Letsgivemore** ♥


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



**KALSTEIN® Climatic Chamber**  
**Environmental Chamber**  
**Model YR05340 // YR05346-1**



**Product Description:**

The Environmental chamber is produced with imported parts with excellent performance, beautiful appearance and good reliability. It is an ideal choice for Laboratory environmental test.

**Performance Characteristics**

- ✓ It has a wide range of temperature and humidity regulation, which can meet various needs of users.
- ✓ It adopts a balanced temperature regulation mode, reach a ideal temperature environment.it has stable balanced heating to carry out high-precision and stable temperature control, timer and over temperature alarm function.
- ✓ The test chamber has the performance of automatically selecting the running cold circuit with the set value of temperature. Simple operation. The equipment has the function of overpressure protection compressor.
- ✓ Big observation window, which is convenient to observe the test state of the sample
- ✓ The double-layer door seal structure ensures that the energy in the cavity is not easy to dissipate. High-low temperature test chamber(YR series): suitable for basic environmental regulations to electrical and electronic products.
- ✓ Test A: Low temperature test method, Test B: High temperature test method. Conduct low temperature and high temperature tests on products. test chamber adopts the most advanced temperature control mode and balanced temperature regulation mode, so it can automatically obtain the reliable and accurate test temperature from high temperature to low temperature or from low temperature to high temperature.
- ✓ Environmental chamber(YR series): It has a wide range of temperature and humidity regulation and control to meet the basic environmental test procedures for electrical and electronic products. Test A: Low temperature test method,
- ✓ Test B: High temperature test method. Test C:humidity test method, Test D:alternative temperature and humidity test method.



**#Let'sgive more** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.




All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



# KALSTEIN® Climatic Chamber

## Environmental Chamber

Model YR05340 // YR05346-1



### Technical Specifications:

Model	Description	Humidity Range	Cooling Range –	Power (W)	Internal Dims (mm)	External Dims (mm)
			Heating up Range			
YR05340	220±10% (-20/0 – 150°C)	20-98%	0.7~1.0°C/min	3500	500x400x400	1420x1000x870
YR05340-1	220±10% (-40 – 150°C)					
YR05342	220±10% (-20/0 – 150°C)	20-98%	–	3500	600x500x500	1500x1050x1130
YR05342-1	220±10%(-40 – 150°C)					
YR05343	220±10% (-20/0 – 150°C)	20-98%	1~3°C/min	4500	750x600x500	1700x1120x1100
YR05343-1	220±10%(-40 – 150°C)					
YR05344	380±10% (-20/0 – 150°C)	20-98%		6500	850x800x600	1860x1320x1270
YR05344-1	380±10% (-40 – 150°C)					
YR05345	380±10% (-20/0 – 150°C)	20-98%		8000	1000x1000x800	1880x1520x1460
YR05345-1	380±10% (-40 – 150°C)					
YR05346	380±10% (-20/0 – 150°C)	20-98%		3000	400x350x350	1420x560x860
YR05346-1	380±10% (-40 – 150°C)					



#### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación  
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through the use of lenses or a system of mirrors that collect and focus light from the object being observed.

**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator is a device that is one of the most common pieces of equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

