




# LABORATORY SECTOR

PATHOLOGICAL ANATOMY



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Pathological Anatomy**  
**Flotation Workstation (Water Bath-Slide Dryer)**  
**Model YR457**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



**KALSTEIN® Pathological Anatomy**  
**Flotation Workstation (Water Bath-Slide Dryer)**

**Model YR457**



**Product Description:**

- ✓ Pure-green digital screen for real-time monitoring of heating temperature, direct and clear, easy to operate; all parameters including preset temperatures, working temperatures, and wording status are real-time displayed
- ✓ Temperature is automatically program-controlled by single-chip microprocessors
- ✓ This device features multiple functions and its easy setting operation can meet the needs of different users
- ✓ DC low-voltage illuminating system and removable transparent heating dish: easy operation and convenient observation
- ✓ New heating mechanism using a new-type high-thermal-conductivity heating element provides even and quick heating: PID-controlled, Triple temperature controls, long lifespan, safe, reliable and energy-saving
- ✓ Temperature is precisely measured by temperature-sensing integrated blocks made in USA, and all settings are automatically stored in the system
- ✓ Special temperature-control and monitoring system characterized by strong resistance to abrasion and corrosion to ensure temperature precise
- ✓ Automatic memory and restoration functions: After startup, all preset temperature data are automatically stored in the system
- ✓ Temperatures for floating and drying slides are individually controlled

**Technical Specifications:**

<b>Model</b>	YR457
<b>Temperature range of water bath</b>	RT – 70°C; automatically maintained
<b>Temperature range of slide dryer</b>	RT – 100°C; automatically maintained
<b>Automatic memory and restoration functions</b>	After startup, all preset temperature
	data are automatically stored in the system.
<b>Temperature Control Precision</b>	±2°C
<b>Dimensions of water bath dish</b>	210×170×60mm(W×D×H)
<b>Dimensions of drying station(50pcs slides)</b>	250×108mm(W×D×H)
<b>Working Voltage</b>	AC 220V±10% 50Hz (standard model)
<b>Power</b>	400W
<b>Dimensions</b>	460×310×135 mm (W×D×H)
<b>Net weight</b>	6kg



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

