




# LABORATORY SECTOR

## BIOSECURITY CABINETS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Biosecurity Cabinets**  
**Horizontal Air Supply Clean Bench**  
**Model YR05733**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



**KALSTEIN® Biosecurity Cabinets**  
**Horizontal Air Supply Clean Bench**  
**Model YR05733**



**Product Description:**

This is a Kind of new type horizontal unidirectional clean bench widely used in electron, national defense, precise instrument and apparatus, pharmacy, chemical experiment and so on

**Main Characteristic:**

- ✓ Adopting arbitrary positioning sliding door system
- ✓ The shell is welded by cold-rolled plates. the surface is electrostatically strayed, and the work surface is SUS304 brushed stainless steel, which is corrosion-resistant and easy to clean
- ✓ The lighting and sterilization systems are interlocked safely, and have the timing function of the ultraviolet lamp
- ✓ Digital display LCD control Interface, more user-friendly design
- ✓ Equipped with a HEPA filter and a primary filter for preliminary filtration, which can effectively extend the service life go the high-efficiency filter

**Technical Specifications:**

<b>Model</b>	YR05733
<b>Air Supply</b>	Vertical Air Supply
<b>Cleanliness level</b>	Class 100@≥0.5μM (US Federal 209E)
<b>Number of colonies</b>	≤0.5pcs/dish·hour (Φ90mm culture plate)
<b>Average Wind Speed</b>	0.25~0.45m/s (Fast and slow)
<b>Noise</b>	≤62dB (A)
<b>Vibration half peak</b>	≤5μM (x、y、z direction)
<b>Illuminance</b>	≥300Lx
<b>Power Supply</b>	AC single-phase 220V/50Hz
<b>Max. Power</b>	800 W
<b>Consumption</b>	
<b>Working area</b>	
<b>W<sub>1</sub>×D<sub>1</sub>×H<sub>1</sub></b>	1 1 40 x 580 x 700
<b>Overall dimension</b>	1 300 x 760 x 1 760
<b>W×D×H</b>	
<b>Specification and quantity of high efficiency filter</b>	
<b>Specification and quantity of fluorescent lamp/ultraviolet lamp</b>	20W x 1 / 36W x 1
<b>Applicable Persons</b>	Double person-single side



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación  
Jacinto Convit**  
<https://www.jacintoconvit.org/>


 **ONE  
TREE  
PLANTED**

<https://onetreepanted.org>

 **Humatem**  
Solidarité et coopération  
biomédicale  
<https://www.humatem.org/>

  
**MANIAPURE  
FOUNDATION USA**  
<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscope is an instrument that uses a series of lenses to magnify images of small objects. This is a common tool in laboratories and research centers.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator is a device that is one of the most common equipment that is used in laboratories. It is a controlled environment (refrigerated space), which allows you to store samples.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

