



# MEDICAL SECTOR

## INFANT RADIANT WARMER



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



# INFANT RADIANT WARMER

## KALSTEIN® Infant Radiant Warmer Infant Radiant Warmer Model YR02190



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Infant Radiant Warmer Infant Radiant Warmer Model YR02190



### Product Description:

Our YR02190 unit offers the best care with the optimal thermoregulation for newborn patients when they are in emergency or therapy jaundice in open care system.

### Features:

- ✓ Pre-warm, manual control and baby skin temperature mode
- ✓ Microprocessor based servo controlled temperature system with far-infrared heater
- ✓ Temperature deviation can be observed and amended in front panel for error measurement
- ✓ Self-checking function for safety affirmative
- ✓ The inclination of bassinet can be adjusted.
- ✓ Heater head can be adjusted horizontally
- ✓ Foldable transparent protector
- ✓ Skin temperature sensor fall protect function to avoid over warming
- ✓ Audible and visible alarm indications for multiple failures with 3 priorities
- ✓ Observe light angle can be adjusted
- ✓ 40% Constant power warming function
- ✓ APAGER timer
- ✓ RS-232 connector

### Application:

- ✓ As our standard designed radiant warmer, the unit can easily regulate the installation and convenient access to the baby patients.
- ✓ Optional accessories
- ✓ Oxygen therapy system unit, oxygen hood
- ✓ Suction unit
- ✓ Resuscitation bag
- ✓ infant head fixating
- ✓ Electrical height adjustment



Model	YR02190
Temperature control mode	Preheating, manual control and skin temperature
Rotary radiant heater	●
Reclining bed	●
Fault alarm	●
APGAR timer	●
X-ray cassette	●
Bed dimensions	big
Superior phototherapy units	○
Lower Light Therapy Units	○
Adjustable working height	N / A
Low pressure suction unit	○
Oxygen therapy unit	○
Oxygen hood	○
Baby headrest	○



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
Microscopes are instruments used for observing microscopic objects. This is a sensitive device that requires a controlled environment to maintain its performance.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigeration is a key element of the most advanced equipment that is used in laboratories. It is a controlled environment (refrigerated space), without fluctuations in temperature.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

