

LABORATORY SECTOR

MOISTURE ANALYZER







All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S.A.S





OUR SERVICES Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💙

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





KALSTEIN® Moisture Analyzer Infrared Moisture Analyzer Model YR05556





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



All rights reserved © KALSTEIN France S. A. S.,

www.kalstein.eu

KALSTEIN® Moisture Analyzer Infrared Moisture Analyzer Model YR05556



Product Description:

- ✓ Infrared lamp/HBM sensor
- ✓ Aluminum shell/Stainless steel chamber
- ✓ High definition LCD
- ✓ Temperature/Time set
- ✓ The percentage of moisture content/dry residual
- ✓ Store historical set
- ✓ Option: Aluminum box/Printer/Grinder

Technical Description:

Model	YR05556
Capacity	110g
Readability	0.001g
Operating Temp	5°C-35°C
Display	LCD HD Back Display
Calibration	External Calibration
Interface	RS232(Standard) / USB(Optional)
Pan Size	Ø90mm
Packing Size	490x350x360mm
Gross Weight	7.5kg
Heating Source	Infrared Lamp
Temp. Sensor	PT100
Temp. Range	40°C-160°C
Temp. Step	1°C
Moisture Range	0.00%-100.00%
Moisture Readability	0.01%
Dry Range	100.00%-0.00%
Dry Readability	0.01%
Temp Setting	40°C-199°C By 1°C Step
Time Setting	1-99min By 1s
Number Of Storage	15
Number Of History	5
Heater	220V ±15% 50HZ/110V±15% 60HZ/400W
Balance	Input:220V 50HZ/110V 60HZ Output:9V



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org





ANY DOUBT? Contact us!



Articles Kalstein (Social Networks and Official Websites)





PARIS - FRANCE **HEADQUARTER** 2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

