




LABORATORY LINES

LAMINAR FLOW
CABINETS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Laboratory Fume Hood, acid and alkali resistance Fume Cupboard Hospital Laboratory equipment



Model Series YR05817 // YR05819

Product Description

- ✓ Front Panel
- ✓ Liner and guide baffles 5mm thick anti melamine board
- ✓ (Window) 5mm thick tempered glass
- ✓ Laboratory dedicated PP Cup
- ✓ Built.in drain pipe
- ✓ Lower Side Plate
- ✓ Lower back plate
- ✓ Lower Movable floor
- ✓ Built in axial fan, air collector and anti-corrosion hose
- ✓ With wind regulator
- ✓ Upper side panel (1.2mm thick galvanized steel plate, epoxy sprayed outside)
- ✓ Built in LED tube
- ✓ Fan and lamp switch
- ✓ Multifunctional splash proof socket
- ✓ 12.7mm corrosion proof solid core, physical and chemical board mesa
- ✓ Door frames and door panels
- ✓ Adjusting foot



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Laboratory Fume Hood, acid and alkali resistance Fume Cupboard Hospital Laboratory equipment



Model Series YR05817 // YR05819



Technical Specifications

Model	YR05817	YR05818	YR05819
External Size (W*D*H) (mm)	1200*800*2350	1500*800*2350	1800*800*2350
Internal Size (W*D*H) (mm)	950*700*1150	1250*700*1150	1550*700*1150
Work Surface Height	800mm		
Max Opening of window	750mm		
Air Velocity/Exhaust air volume	0.4-0.6m/s,800-1000m³/h		
Noise	≤60dB		
LED Lamp	10W*1		
Blower	Built-in axial flow blower;Speed adjustable.		
Front Window	Manual,5mm toughened glass,height adjustable.		
Power Supply	AC220V±10%,50/60Hz;		
Consumption	100W		
Material	Main Body	1.2mm thickness cold-rolled steel with epoxy resin powder coating.	
	Inner Side	5mm thickness compact,resistant to strong acid and alkali.	
	Work Table	Chemical resistant phenolic resin, thickness 12.7mm.	
Standard Accessory	LED lamp,Built-in axial flow blower with speed adjustable.		
	Waterproof socket*3		
	3 meter PP duct,Diameter:250mm		
Optional Accessory	epoxy resin board or ceramic board.,Water tap,Gas tap,Water Sink		
	Outer centrifugal blower(300W),resistant to strong acid,alkali.		
Gross Weight	235	265	300



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through the use of lenses or mirrors that collect and focus light from the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most common pieces of equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog


Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

