




LABORATORY LINES

OIL BATHS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Water Baths

Oil Bath

Series YR04973

YR oil bath series is widely used in bioengineering, medicine, food, chemical industry, gold treatment, petroleum and other fields. To provide users with high-precision, controlled, temperature uniform constant place. Also used in research institute, colleges and universities, factory laboratory, inspection department, etc.

Features:

- ✓ Microcomputer intelligent control system has the advantages of rapid temperature rise, stable temperature and easy operation
- ✓ Dual purpose water and oil: the maximum temperature can reach 250°C
- ✓ Led dual window digital display temperature measurement value and temperature setting value respectively, convenient to touch the key
- ✓ It is suitable for chemical reaction, biomedicine and other industries which need high temperature field
- ✓ Optional magnetic stirring device for high temperature reaction



Model (Integrated chamber)	YR04973	YR04974	YR04975	YR04976
Controller			PID	
Voltage			220V / 50Hz 110V/60Hz	
Power (W)	1 000W	1 200W	1 600Ww	2 000W
Temperature – Range (°C)			RT. +5~250°C	
Temperature Resolution			0.1 °C	
Temperature Fluctuation			≤±0.5°C	
Division Accuracy			≤±1°C	
Timer			999 minutes	
Recycle Design			Optional	
Magnetic			Optional	
Safety	Over current protection, Over temperature alarm; Mechanical temperature limiter			
Chamber Size (mm)	300 x 240 x 150	300 x 240 x 200	500 x 290 x 150	500 x 290 x 200
Product Size (mm)	390 x 280 x 240	390 x 280 x 290	600 x 350 x 240	600 x 350 x 290
Packing Size (mm)	440 x 330 x 290	440 x 330 x 340	650 x 400 x 290	650 x 400 x 340



#Letsgivemore ❤️

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

