




LABORATORY SECTOR

PH METERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® PH Meters
Pocket TDS Meter
Model YR01813



#Letsgivemore ♥


Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® PH Meters
Pocket TDS Meter
Model YR01813



Product Description:

Economical pocket TDS tester, 1 point calibration. The meter is suitable for measuring the low concentration liquids, accuracy: 1% F.S.

Features:

- ✓ Low-range pocket TDS tester allows adjusting the TDS conversion factor.
- ✓ Platinum conductivity cell provides the quick and reliable measurement results.
- ✓ 1 point push-button calibration ensures accurate measurement over the entire range.
- ✓ Hold function momentarily freezes reading for easy viewing and recording.
- ✓ Auto-Power Off effectively conserves battery life.
- ✓ Replaceable electrode module reduces maintenance and replacement cost.
- ✓ Waterproof feature to ensure complete protection in harsh environment.

Applications:

- ✓ Water and wastewater treatment
- ✓ Aquaculture and other industrial applications



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® PH Meters

Pocket TDS Meter

Model YR01812



Technical Description:

Model	YR01813
TDS	
Range	0.05 to 10.00ppt
Resolution	0.01ppt
Accuracy	±1% F.S.
Calibration Point	1 point
Calibration Solution	7447ppm
Temperature	
Range	0 to 50°C
Resolution	1°C
Accuracy	±1°C
General Specifications	
Temperature Compensation	0 to 50°C, automatic
TDS Factor	0.4 to 1.0 (default 0.5)
Hold Function	Manual
Auto-Power Off	8 minutes after last key pressed
Sensor Type	ECSCAN-C1-100K
Display	Single-line LCD (21×21mm)
Power	3×1.5V LR44 micro alkaline batteries
Battery Life	Approximately 150 hours of continuous use
Dimensions	185(L)×40(Dia.)mm
Weight	100g



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

