



LABORATORY SECTOR

BALANCES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Balances
Semi-Micro Analytical Balance
Model YR05562 // YR05567



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Balances
Semi-Micro Analytical Balance
Model YR05562 // YR05567



Product Description:

- ✓ Recipe database: 99 data's can be stored
- ✓ Automatic internal calibration
- ✓ Tare/counting/percentage test/animal weighing
- ✓ Baud rate is adjustable/ Density function
- ✓ Unit conversion (g/mg/ct/oz)...
- ✓ Multi-Language : English/
German/French/Italian/Spanish/Portuguese
- ✓ Keep weight values in GLP/ISO format, Includes
date, time, and serial number
- ✓ Plastic protection cover for keypad
- ✓ Sensitivity/speed adjustable
- ✓ Overload alarm/fault alarm/level indicator

Optional Accessories:

- ✓ ION-A15 Static Eliminator
- ✓ Balance plate system. USB version (BL0329) or
Bluetooth version (BL0330)
- ✓ DPP-250-BT Bluetooth printer for system (AC013)
- ✓ TLP-50 serial printer with date/time (C054)
- ✓ DPP-250 serial printer (AC007)
- ✓ STAT serial printer with statistical function (AC022)
- ✓ Alphanumeric external keyboard (AC005)
- ✓ Serial to USB converter (E1002)
- ✓ Serial cable for serial output to printer or computer
(E743)
- ✓ Factory calibration certificate (BL0333)



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Balances
Semi-Micro Analytical Balance
Model YR05562 // YR05567



Technical Specifications:

Model	YR05562	YR05563	YR05564	YR05565	YR05566	YR05567
Capacity(g)	62g/120g	82g/220g	102g/220g	122g	162g	210g
Resolution(g)	0.01mg/0.1mg			0.01mg		
Min weighing(g)	60mg		100mg	100mg		
Stable Time	≤6s					
Display	LCD (white back light with black font)					
Pan size	Φ80mm					
Operation temp	5-25°C 20-80% humidity					
Repeat ability	0.03mg/0.1mg	0.05mg/0.1mg		0.05mg	0.05mg	
Liner	±0.08/±0.2mg	±0.1/ ±0.2mg		±0.1mg	±0.1mg	
Draft shield size	152mm×160mm×225mm					
Cal.weight	Internal Calibration					
Option	RS232C & RS485、Printer					
Dimension	345×215×345mm					
G.N . weight	8.25kg/5.8kg				8.3kg/5.85kg	
IONISER	Option				Standard	
Language	Italian/English/Portuguese/Deutsch/France/Spanish					
Packing Size	510x400x620mm(No pallet)					
Power Supply	110-230V 50/60Hz Out DC 12V 1A					



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>

 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscopes in a laboratory are one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

