



LABORATORY SECTOR

REFRIGERATORS AND FREEZERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Refrigerators and Freezers
Sprayed Plate Steel Medical Chest Freezer

Model YR05318



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers

Sprayed Plate Steel Medical Chest Freezer

Model YR05318



Safety

- ✓ Build-in back-up battery to display temperature and to operate audible and visual alarm system for up to 12 hours without AC power.
- ✓ Visual and audible alarm system: High/low temperature, Sensor error, Thermostat failure alarm, power failure alarm
- ✓ Power failure protection: turn on delay of the cooling system when power failure, restart delay protection.
- ✓ Keyboard lock, password protection to avoid set parameter randomly.
- ✓ Wide available voltage range: 198V~242V.

Special:

- ✓ 4 visual and audible alarms. 1 Standard test hole(25mm).



Product Description:

485L Chest Medical Freezer From China Made Of Sprayed Plate Steel For Hospital And Laboratory Use

Promed Medical Freezers are Designed to Store Bio-Samples Safety

- ✓ Build-in back-up battery to display temperature and to operate audible and visual alarm system for up to 12 hours without AC power.
- ✓ Visual and audible alarm system: High/low temperature, Sensor error, Thermostat failure alarm, power failure alarm
- ✓ Power failure protection: turn on delay of the cooling system when power failure, restart delay protection.
- ✓ Keyboard lock, password protection to avoid set parameter randomly.
- ✓ Wide available voltage range: 198V~242V.

Humanized Design:

- ✓ 4 casters,2 casters with stopper for easy moving and fixation.
- ✓ 1 Standard test hole(25mm),convenient for temperature monitoring.

Features:

- ✓ Wide temperature range :10°C~32°Cambient temperature can reach -10°C~-25°C
- ✓ Accurate Temperature control technology to make the temperature maintain-10°C~-25°Cand accuracy 0.1 °C.
- ✓ SPCC sprayed steel external material and inner chamber for friendly using and easy cleaning.
- ✓ High-density polyurethane foam.

Certificate:

- ✓ ISO14001.ISO9001.ISO13485.ISO45001
- ✓ FDA List and CE

Warranty:

- ✓ 1 year free spare parts and remote technical support.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers

Sprayed Plate Steel Medical Chest Freezer

Model YR05318



Technical Description:

Model	YR05318
Performance	
Temperature Range	-10~-25°C
Ambient Temperature	+10~+32°C
Controller	Microprocessor
Display	Digital Display
Refrigeration System	
Compressor/Brand	1 /CUBIGEL
Refrigerant	R290
Cooling Type	Direct cooling
Defrost Type	Manual
Electrical Data	
Power Supply(V/Hz)	220~240V/50~60;115/60
Rated Power(W)	295
Electrical Current(A)	2.45
Power Consumption	3.45 kWh/24h
Dimensions	
Capacity	485L/17.12 cu.ft
Interior Size(W•D•H)	1455*525*665mm
Exterior Size(W•D•H)	1655*755*915mm
Package Size(W•D•H)	1710*885*1020mm
N.t./G.t.(approx)	124/138KG
Loading Qty(20"/40"/40"H)	12/26/26
Alarm	
High/Low Temperature	Y
Sensor Error	Y
Power Failure	Y
Accessories	
Caster	6
Remote alarm port	/
Test hole	1/25MM



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

