




MEDICAL SECTOR

INFUSION PUMP



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Infusion Pump
Ultrasound Laptop 15 Inch Touch Screen
Model YR05165



#Letsgivingore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Infusion Pump Ultrasound Laptop 15 Inch Touch Screen Model YR05165



Product Description:

- ✓ High-resolution 15.1" color LED backlight display with high contrast and wide viewing angle, touch screen is optional
- ✓ Display mode: B, 2B, 4B, B / M, M
- ✓ Image storage: 4G hard drive to permanently store around 5000 image frames
- ✓ One Key Store, One Key Review, One Key Print
- ✓ One-key transfer of current frozen image to network workstation to complete your image text report and print directly
- ✓ The built-in battery can work about 3 hours
- ✓ software to test General, Obstetrics, Gynecology, Vascular, Cardiac, Small Organ, MSK, Joint, etc.
- ✓ Multilanguage function: Spanish, English, French, Russian and Portuguese



Technical Description:

Model	YR05165
Menu operation system, with different languages.	according to requirement including Spanish, English, French, Russian and Portuguese
Body Brands	total of 64 kinds of body markings indicating the position of the probe
Display Mode	B, 2B, 4B, B / M, M;
Adjustable Depth	10 degree, max depth 236mm
To win	0 ~ 127dB
Contrast	27 ~ 90dB
Grayscale	256, 8-step smart TGC
Frame Rate	30 cond boxes, 8 kinds of pseudo colors
Image Storage	4G hard drive to permanently store around 5000
Marcos images	no loss when power off. 256 frames; Puncture guide function with 2
Cinema Loop	correctable puncture guide lines (angle and position adjustable); the gravel software package with real-time position line measurement
Normal Measure	Distance, perimeter, area, volume angle, length and area stenosis ratio measurement.
Cardiovascular Measurement	Depth, slope, heart rate, cycle



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Infusion Pump
Ultrasound Laptop 15 Inch Touch Screen
Model YR05165



Cardiovascular Measurement	Depth, slope, heart rate, cycle
gynecology measurement	and analysis:
	general gynecological measurement and gynecological report
Obstetric Measurement	including: 15 kinds of obstetric software,
	various GA measurements,
	amniotic fluid index,
	fetal growth curve and OB report.
	residual urine volume,
Urology Measurement and Analysis	prostate volume and
	urology report.
	concise and practical
One Key Store	the current frozen image to the network workstation to
One Key Transfer	complete your image text report and print directly;
	Dual TV output: PAL / NTSC, you can connect a video printer when necessary;
Adjustable Display Brightness	
Power supply	100-240V ~ 1.2-0.6A frequency: 50-60Hz;
Power Supply Adapter Output	DC 12.8 V 3.0 A
Main Unit Weight	2.9 kg 9 (without accessories)
Main Unit Dimension	380 × 300 × 65mm



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

