




Kalstein®

A different accompaniment, at your service

LABORATORY SECTOR

MICROSCOPE



All rights reserved © KALSTEIN France S. A. S.,



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Microscope

Zoom Stereo Microscope, 0.67-4.5x,1:6.7

Model YR0A23.3667L



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



YR0A23.3667L Details



Eyeiece WF10X/20mm, High Eyepoint, Eyeiece Tube Diopter Adjustable



0.67x-4.5x, Magnification Up To 3.35x-270x With Optional Eyeieces And Auxiliary Objectives



YR0A54.3640 35 Kinds Stands Optional For Different Usage & Condition



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



YR0A23.3667L Specification



YR0A23.3667



YR0A23.3667L



YR0A23.3667L Series Zoom Stereo Microscope Specification


		B	T
Head	Binocular Head, 45° Inclined, 360° Rotatable, Interpupillary Adjustment 54-75mm	•	
	Trinocular Head, 45° Inclined, 360° Rotatable, Interpupillary Adjustment 54-75mm, With Eyepiece Adapter Focus Adjustable, Light Split Fixed P50:E50, Simul Focal Design To Show Image On Both Eyepiece & Photo Port At Same Time		•
Eyepiece	WF10x/22mm, Dia.30mm, High Eyepoint, Diopter Adjustable		
Zoom Range	0.67x~4.5x		
Zoom Ratio	1:6.7		
Magnification	6.7~45x, Up To 2x-270x With Optional Eyepieces & Auxiliary Objectives		
Working Distance	105mm, Up To 26~287mm With Auxiliary Objectives		
Field Of View	0.9-109.5mm (Detail See Optical Data Sheet)		
Focusing	Pole Stand Coarse Focusing Range 50mm, Tension Adjustable Track Stand Coarse Focusing Range 74mm, Tension Adjustable		
Stand	35 Kinds Stand Optional, Details See A54.3630 Choose Model According To Cata. No. Like: YR0A23.3667L-B1 Means Binocular Head With B1 Stand YR0A23.3667L-B1T Means Trinocular Head With B1 Stand		

YR0A23.3667L Series Zoom Stereo Microscope Optical Data Sheet

Zoom Lens		0.67x~4.5x						
Auxiliary Lens		0.3x	0.5x	0.7x	0.75x	1.0x(Standard)	1.5x	2.0x
Working Distance		287mm	177mm	120mm	117mm	105mm	47mm	26mm
WF10x/22mm	Magnification	2.0-13.5x	3.4-22.5x	4.7x-31.5x	5.0-33.8x	6.7-45.0x	10.1-67.5x	13.4-90.0x
	View Field(mm)	109.5-16.3	65.7-9.8	46.9-6.9	43.8-6.5	32.8-4.9	21.9-3.3	16.4-2.4
WF15x/16mm	Magnification	3.0-20.3x	5.0-33.8x	7x-47.2x	7.5-50.6x	10.1-67.5x	15.1-101.3x	20.1-135x
	View Field(mm)	79.6-11.9	47.8-7.1	34.1-5.1	31.8-4.7	23.9-3.6	15.9-2.4	11.9-1.8
WF20x/12mm	Magnification	4.0-27.0x	6.7-45.0x	9.3x-63x	10.1-67.5x	13.4-90.0x	20.1-135.0x	26.8-180x
	View Field(mm)	59.7-8.9	35.8-5.3	25.6-3.8	23.9-3.6	17.9-2.7	11.9-1.8	9-1.3
WF25x/9mm	Magnification	5.0-33.8x	8.4-56.3x	11.7x-78.7x	12.6-84.4x	16.8-112.5x	25.1-168.8x	33.5-225x
	View Field(mm)	44.8-6.7	26.9-4.0	19.2-2.8	17.9-2.7	13.4-2.0	9.0-1.3	6.7-1.0
WF30x/8mm	Magnification	6.0-40.5x	10.2-67.5x	14.1x-94.5x	15.0-101.4x	20.1-135x	30.3-202.5x	40.2-270x
	View Field(mm)	37.3-5.6	23.8-3.6	16-2.25	14.1-2.25	11.9-1.8	7.9-1.2	6.0-0.9

Optional Accessories

Head	Binocular Head + WF10x/20mm	
	Trinocular Head + WF10x/20mm	
Eyepiece	WF10x/22mm, High Eyepoint, Diopter Adjustable, Dia.30mm	
	WF15x/16mm, High Eyepoint, Diopter Adjustable, Dia.30mm	
	WF20x/12mm, High Eyepoint, Diopter Adjustable, Dia.30mm	
	WF25x/9mm, High Eyepoint, Diopter Adjustable, Dia.30mm	
	WF30x/8mm, High Eyepoint, Diopter Adjustable, Dia.30mm	
Auxiliary Lens	0.3x, W.D. 287mm	
	0.5x, W.D. 177mm	
	0.7x, W.D. 113mm	
	0.75x, W.D. 116mm	
	1.5x, W.D. 47mm	
	2.0x, W.D. 26mm	
CCD Adapter	Eyepiece Adapter Focus Adjustable	
	1.0x, Focus Adjustable	
	0.50x, Focus Adjustable, For 1/2" CCD	
	0.35x, Focus Adjustable, For 1/3" CCD	



YR0A54.3630 35 Kinds Stand Optional

1

YR0A54.3630 35 Kinds Stand Optional, Details See Choose Model According to Catalogue No. Like: YR0A23.3667L-B1
Means Binocular Head With B1 Stand YR0A23.3667L-B1T Means Trinocular Head With B1 Stand



2



A1 Holder



B1 Pole Stand



B1 Stand

3



A3 Holder



STL4 Boom Stand



STL4 Stand



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





Binocular Head

35 Stands



YR0A23.3667L-B1



YR0A23.3667L-B2



YR0A23.3667L-B2L



YR0A23.3667L-B3



YR0A23.3667L-B3LS



YR0A23.3667L-B2LD



#Let'sgive**more**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



MICROSCOPE



YR0A23.3667L-B4



YR0A23.3667L-B5



YR0A23.3667L-B6



YR0A23.3667L-B7



YR0A23.3667L-B8L



YR0A23.3667L-B8LD



YR0A23.3667L-B9L



YR0A23.3667L-B9LD



YR0A23.3667L-B8LS



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



MICROSCOPE

CE



YR0A23.3667L-L



YR0A23.3667L-2L



YR0A23.3667L-B9LS



YR0A23.3667L-J1



YR0A23.3667L-J2



YR0A23.3667L-J3L



YR0A23.3667L-J4L



YR0A23.3667L-J4LD



YR0A23.3667L-J3LD



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



MICROSCOPE



YR0A23.3667L-STL1



YR0A23.3667L-STL2



YR0A23.3667L-STL3



YR0A23.3667L-STL4



YR0A23.3667L-STL5



YR0A23.3667L-STL8



YR0A23.3667L-STL6A



YR0A23.3667L-STL6B



YR0A23.3667L-STL6C



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



YR0A23.3667L Series Trinocular Models



Trinocular Head



35 Stands



YR0A23.3667L-B1T



YR0A23.3667L-B2T



YR0A23.3667L-B2LT



YR0A23.3667L-B3T



YR0A23.3667L-B3LST



YR0A23.3667L-B2LDT



#Letsgivemore Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





YR0A23.3667L-B4T



YR0A23.3667L-B5T



YR0A23.3667L-B6T



YR0A23.3667L-B7T



YR0A23.3667L-B8LT



YR0A23.3667L-B8LDT



YR0A23.3667L-B9LT



YR0A23.3667L-B9LDT



YR0A23.3667L-B8LST



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



MICROSCOPE



YR0A23.3667L-LT



YR0A23.3667L-2LT



YR0A23.3667L-B9LST



YR0A23.3667L-J1T



YR0A23.3667L-J2T



YR0A23.3667L-J3LT



A23.3667L-J4LT



A23.3667L-J4LDT



A23.3667L-J3LDT



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



MICROSCOPE



YR0A23.3667L-STL1T



YR0A23.3667L-STL2T



YR0A23.3667L-STL3T



YR0A23.3667L-STL4T



YR0A23.3667L-STL5T



YR0A23.3667L-STL8T



YR0A23.3667L-STL6AT



YR0A23.3667L-STL6BT



YR0A23.3667L-STL6CT



#Let'sgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





YR0A51.3622-1022
Eyepiece WF10x/22mm



YR0A51.3622-1516
Eyepiece WF15x/16mm



YR0A51.3622-2012
Eyepiece WF20x/12mm



YR0A51.3622-2509
Eyepiece WF25x/9mm



YR0A51.3622-3008
Eyepiece WF30x/8mm



YR0A52.3622-03
0.3x, W.D. 287mm



YR0A52.3622-05
0.5x, W.D. 177mm



YR0A52.3622-07
0.7x, W.D. 113mm



A52.3622-75 0.75x,
W.D.YR0. 116mm



YR0A52.3622-15
1.5x, W.D. 47mm



YR0A52.3622-20
2.0x, W.D. 26mm



YR0A55.3623-E
Eyepiece Adapter



YR0A55.3623-10
1.0x



YR0A55.3623-05
0.50x, For 1/2" CCD



YR0A55.3623-35
0.35x, For 1/3" CCD



YR0A56.2103-A
LED Ring Light, 56 LEDs



YR0A56.1211-B
LED Ring Light



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



YR0A22.36, YR0A23.36 Total 8 Stereo Microscope Series



YR0A22.3660
1x/3x, 2x/4x



YR0A22.3660N
1x/3x, 2x/4x



YR0A22.3670
1x/3x, 2x/4x



YR0A23.3645
0.7x~4.5x, 1:6.5
Binocular, Trinocular



YR0A23.3645
0.7x~4.5x, 1:6.5
Binocular, Trinocular



YR0A32.3645
0.7x~4.5x, 1:6.5
Digital Head



YR0A23.3667
0.67x~4.5x, 1:6.7
Binocular, Trinocular



YR0A23.3667L
0.67x~4.5x, 1:6.7
Binocular, Trinocular



Kalstein[®]

A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>


Articles Kalstein
(Social Networks and Official Websites)



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago

